

# Glycopeptide profiling in a trapped ion mobility quadrupole time-of-flight mass spectrometer (timsTOF Pro)

Wessels HJCT<sup>1</sup>, Bärenfänger M<sup>1</sup>, Post M<sup>1</sup>, Zijlstra F<sup>1</sup>, Marx K<sup>3</sup>, Wei W<sup>4</sup>, Li H<sup>4</sup>, Zhang T<sup>4</sup>, Shan P<sup>4</sup>, Schmit P<sup>2</sup>, van Gool AJ<sup>1</sup>, Lefeber DJ<sup>1</sup>.


<sup>1</sup> Translational Metabolic Laboratory, Radboudumc, Nijmegen, Netherlands

<sup>2</sup> Bruker, Wissembourg, France

<sup>3</sup> Bruker Daltonik GmbH, Bremen, Germany


<sup>4</sup> Bioinformatics Solutions, Waterloo, Canada

# Introduction




Recent introduction of parallel accumulation serial fragmentation (PASEF) on a trapped ion mobility quadrupole time-of-flight mass spectrometer provides unique possibilities for comprehensive glycopeptides profiling in complex samples such as blood plasma. Glycopeptide profiling of patient samples holds great potential to diagnose, monitor and understand both genetic and acquired human diseases. Glycosylated peptides present analytical challenges for LC-MS/MS analysis with respect to chromatographic separation, electrospray ionization, and structural elucidation by collision induced dissociation experiments. In this work we set out to optimize the data generation workflow on the timsTOF Pro platform for N-glycopeptide profiling. Parameter optimizations discussed in this communication are: C18RP mobile phase composition (**slide 3**), TIMS and QqTOF ion optics (**slide 4**), ESI conditions (**slide 5**) and PASEF parameters (**slides 6-8**).

# Improved C18RP separation of glycopeptides by TFA ion pairing



Base peak chromatograms of tryptic digest of model glycoproteins (TRFE\_HUMAN, A1AG1\_HUMAN, FETUA\_BOVIN, Trastuzumab) using 0.1% FA and 0.1% FA, 0.01% additives in H<sub>2</sub>O (solvent A) and ACN (solvent B). Gradient: linear increase of 3-45% solvent B in 15 minutes at 500nl/min. Trapping column: Acclaim C18RP RLSC (0.075 x 20mm, 3µm particles, 100Å pore size; Thermo scientific). Analytical column: C18RP nanoElute FIFTEEN analytical column (0.075 x 150mm, 1.9µm particles, 120Å pore size; Bruker Daltonics) operated at 45°C.




Glycoform	APEX Rt (min)		Full width (sec)		Tailing	
	FA	FA TFA	FA	FA TFA	FA	FA TFA
H5N4S2	3.3	3.46	7.3	5.8	1.403	1.086
H6N5S3	3.61	3.43	14.8	6.7	0.966	1.027
H7N6S4	5.56	3.38	165	6	N.D.	1.282


Extracted ion chromatograms of disialo-, trisialo-, and tetrasialo glycoforms of peptide ENGTISR from human Alpha-1-acid glycoprotein 1 (A1AG1\_HUMAN). With 0.1% FA solvent composition glycoforms are separated according to their respective number of sialic acids. The EIC of the tetrasialo glycoform of ENGTISR shows extreme peak broadening and tailing with 0.1% FA. Addition of 0.01% TFA significantly improves chromatographic separation of all glycoforms but achieves a remarkable 27.5x reduction in full peak width for the tetrasialo glycoform.

Reversed phase C18 liquid chromatography separation of glycopeptides with typical 0.1% formic acid (FA) solvent composition performs poorly for highly sialylated glycopeptides. Glycoforms of the same peptide show increased retention, peak broadening and increased tailing when the number of sialic acids in the glycan structure increases. Addition of low levels of trifluoroacetic acid (TFA) as ion pair reagent mitigates interactions between sialic acids and the stationary phase leading to near co-elution of all glycoforms of a peptide with sharp chromatographic peaks with minimal tailing. Signal intensity due to ion suppression by TFA is only ~50% thanks to acetonitrile nanoBooster dopant.

# TIMS delta potentials and ion optics tuning




TIMS		MS1 and low CE MS2		high CE MS2	
D1	-20V	Transfer time	100 $\mu$ s	Transfer time	64 $\mu$ s
D2	-150V	Collision RF	2000Vpp	Collision RF	1600Vpp
D3	110V	Pre Pulse Storage	13 $\mu$ s	Pre Pulse Storage	10 $\mu$ s
D4	110V				
D5	0V				
D6	55V				




TIMS delta potentials and QqTOF ion optics were tuned for optimal transmission of intact glycopeptides using tunemix and GluFib direct infusion as well as LC-MS analysis of glycopeptides from glycoprotein standards. First we optimized transmission by ramping individual TIMS delta potentials whilst monitoring intact analyte ion transmission as well as detectable fragment ions. GluFib results show that delta 3 and delta 6 TIMS potentials lead to undesired in-TIMS fragmentation when set to high values. Glycopeptides measurements indicate that delta 6 TIMS potential is primarily responsible for in-TIMS fragmentation. QqTOF ion optics were optimized using tunemix for optimal mass range transmission of both intact glycopeptides in MS1 or glycan-moiety fragments in low collision energy MS2 experiments. Finally, we also optimized ion optics parameters for transmission of peptide-moiety fragments generated in high collision energy MS2 experiments. Optimal key parameters are listed in the table.

# nanoBooster ESI nitrogen source gas dopants



CaptiveSprayer with nanoBooster image and bar graph of summed peak areas of glycopeptides from a tryptic digest mixture of TRFE\_HUMAN, A1AG1\_HUMAN, FETUA\_BOVIN and Trastuzumab using different nanoBooster dopants. Acetonitrile provides highest possible glycopeptide signal intensities and highest precursor ion charge states for optimal CID fragmentation of peptide-moieties.




Extracted ion mobilograms for selected glycoforms from the tetrasialo N-glycopeptide ENGTISR from A1AG\_HUMAN using different nanoBooster dopants.


Glycopeptides typically suffer from poor ionization efficiency in electrospray ionization. Use of dopant enriched N<sub>2</sub> gas in the captive sprayer ESI source significantly enhances the ionization efficiency of glycopeptides but also results in supercharging (acetonitrile) or subcharging (primary alcohols) of ions depending on the solvent used. NanoBooster solvents have a significant impact on TIMS separation as the increased gas density leads to an increased collisional force that pushes ions deeper into the TIMS funnels. Positioning of ions within the TIMS funnels was corrected by reducing the gas velocity and polyalanine as collisional cross section standards was used to calibrate reduced mobility values (1/KO). In general, only minor shifts in 1/KO values could be observed for the vast majority of glycopeptides with no notable differences in mobility peak shape, number of mobility features or resolution. However, sporadic cases could be observed where both 1/KO as well as the mobility peak shapes of glycopeptide ions were significantly influenced in a dopant specific manner.

# Optimization of PASEF parameters


### Detected 4D features



### Acquired MS2 spectra



### Unique glycopeptide identifications







Results obtained for enriched glycopeptides from control plasma at various accumulation + ramp times in combination with various target PASEF intensities. Glycopeptides were separated using 60 minute linear gradients of 3-45% acetonitrile in 0.1% FA, 0.01% TFA at 500nl/min. 4D feature detection was performed using Peaks studio X+ (BSI). MGF files were generated using DataAnalysis 5.2 and analyzed using BioPharma Compass 3 (Bruker Daltonics). MS2 spectra were classified by ProteinScape using a HexNAc-HexNAc-Hex-Hex N-glycan core fragmentation pattern and subsequently subjected to glycan database searches using GlycoQuest (Bruker Daltonics) and protein database searches using MASCOT 2.5.1 (Matrix Science).

Using optimized low- and high-collision energies for respective glycan- and peptide-moiety fragmentation we performed measurements with combinations of increasing accumulation + ramp times with different MS2 target intensities for PASEF. The number of detected MS1 features that serve as precursors for PASEF steadily decreases with lower accumulation and ramp times due to the abundance distribution of glycoproteins in plasma. As a result the total number of acquired MS2 scales accordingly with increased accumulation + ramp times as more precursor are detected. The target PASEF intensity scales the number of replicate MS2 spectra recorded and summed for a single precursor according to its precursor intensity in MS1. Higher target PASEF intensities increase MS2 sensitivity but negatively impact the number of unique precursors that can be fragmented per duty cycle. Best results for plasma glycoproteomics were obtained by setting parameters that maximize the number of ions accumulated in the TIMS section and aim for high MS2 sensitivity.

# Glycopeptide ID example: dominant glycoforms of hemopexin

10	20	30	40	50	60
MARVLGAPVA	LGLWSLCWSL	AIATPLPPTS	AHGNVAEGET	KPPDPVTERC	SDGWSFDATT
70	80	90	100	110	120
LDDNGTMLFF	KGEFVWVSKH	WDRELISERW	KNFPPSPVDA	FRQGHNSVFL	IKGDKVWVYP
130	140	150	160	170	180
PEKKEKGYPK	LLQDEFPGIP	SPLDAAVECH	RGECQAEGLV	FFQGDREWFV	DLATGTMKER
190	200	210	220	230	240
SWPAVGNCS	ALRWLGRYYC	FQGNQFLRFD	PVRGEVPPRY	PRDVRDYFMP	CPGRGHGHRN
250	260	270	280	290	300
GTGHGNSIHH	GPEYMRCSPH	LVLSTALSDN	HGATYAFSGT	HYWRLDTSRD	GWHSWPIAHQ
310	320	330	340	350	360
WPQGSAPVDA	AFSWEELVYL	VQGTQVYVFL	TKGGYTLVSG	YPKRLEKEVG	TPHGIIIDSV
370	380	390	400	410	420
DAAFICPGSS	RLHIMAGRRL	WLDLKSQAQ	ATWTELPWPH	EKVDGALCME	KSLGPNSCSA
430	440	450	460	470	
NGPGLYLIGH	PNLYCYSDE	KLNAKALPQ	PQNVTSLLGC	TH	




# Peaks GlycoFinder: prototype 4D glycoproteomics software

PEAKS GlycoFinder is currently being developed by Bioinformatics Solutions Inc. In collaboration with Bruker Daltonics and Radboudumc for comprehensive ion mobility aware glycoproteomics. The current prototype combines novel data processing with advanced glycan- and peptide-moiety identification strategies using both database searches and *de novo* sequencing for deep glycoproteome profiling.

4D feature detection for label-free quantitation combined with bioinformatic data processing enables users to explore protein glycobiochemistry through intuitive graphical representations of the data from different perspectives.

Preliminary results for a small subset of 6 full plasma glycopeptide data files acquired on a timsTOF Pro with optimized acquisition parameters yielded 2423 unique glycopeptides from 187 glycoproteins. Further optimization of algorithms is in progress to further improve glycopeptide identification results.



Screenshot of the prototype PEAKS GlycoFinder user interface




# Conclusions

We have characterized and optimized critical parameters for successful implementation of N-glycopeptide profiling in human plasma on the timsTOF Pro platform. Follow up work is now focussed on glycopeptide software development and application of PASEF glycopeptide technology in clinical applications to improve patient care.

## Take home messages:

- Addition of TFA to mobile phase is required for analysis of highly sialylated glycoforms using ReproSil C18AQ2 stationary phase
- High TIMS delta potentials lead to in-TIMS fragmentation of glycopeptides
- nanoBooster N2 source gas dopants are compatible with TIMS, significantly increase sensitivity, and shift charge state distributions
- Best plasma N-glycopeptide identification results are obtained with PASEF parameters that take advantage of the timsTOF Pro MS/MS acquisition speed to maximize sensitivity
- Application of our optimized PASEF workflow in combination with prototype 4D glycoproteomics software enables the analysis of thousands of N-glycopeptides in unfractionated human plasma samples



# Acknowledgements

**Radboudumc**  
university medical center

## Translational Metabolic Laboratory

Alain van Gool  
Anouk Suppers  
Dirk Lefeber  
Esther Willems  
Fokje Zijlstra  
Hans Wessels  
Jolein Gloerich  
Koen Rademaker  
Maurice van Dael  
Melissa Bärenfänger  
Merel Post  
Monique van Scherpenzeel  
Nurulamin Bin Abu Bakar  
Purva Kulkarni  
Ron Wevers  
Thijs Deenik  
Yanniek Lagarde



## Bruker Daltonics

Gary Kruppa  
Kristina Marx  
Markus Lubeck  
Pierre-Olivier Schmit  
Scarlet Beck  
Stuart Pengeley



## Bioinformatics Solutions inc.

Hui Li  
Paul Shan  
Tony Zhang  
Winter Wei