PREOM/CS

# Overcoming the dynamic range difficulties in CSF samples through a novel enrichment step for biomarker discovery studies

Katharina Limm<sup>1</sup>, **Zuzana Demianova<sup>1</sup>**, Katrin Hartinger<sup>1</sup>, Nils A. Kulak<sup>1</sup> <sup>1</sup>PreOmics GmbH, Martinsried, Germany

## SPOTLIGHT

- simple-to-use and optimized • Fast, protocol for the preparation of complex samples like cerebrospinal fluid, in proteomic profiling high-throughput approach for biomarker discovery
- ENRICH-iST tackles the challenge of a high diluted protein content and a high dynamic abundance range, which is typical for CSF
- Streamlined solution to reduce the high dynamic range for greater proteomic depth
- Deeper with coverage proteome biological information for preserved additional potential discovering biomarkers



## **MATERIALS & METHODS**

#### Source:

Cerebrospinal fluid samples from rat.

### Sample preparation:

10 µL of CSF followed by iST-BCT protocol (PreOmics) for neat CSF and 20-200 µL of CSF followed by ENRICH-iST protocol.

### LC-MS analysis:

NanoELUTE2 (24 SPD method) flow rate of 300 nl/min on a reversed-phase C18, IonOpticks Aurora column (25 cm x 75  $\mu$ m, 1.7  $\mu$ m) coupled to a timsTOF HT in dia-PASEF mode. Sample load was adjusted to 300 ng load on column.

### Data processing:

Burker ProteoScape<sup>™</sup> was used for real-time processing of each single MS-run. Afterwards post processing was performed assigning single runs to sub-goups of replicates without match -between-runs.

groups could be obtained against neat CSF.





Additinal identified proteins in ENRICH-iST prepared samples. Comparing the IDs found in differently regulated proteins detected in both conditions (neat and the neat and ENRICH-iST workflows, almost all proteins were detected by ENRICH-iST. Further ENRICH-iST). In both comparisons, proteins with high abundance were not identifications were possible with a higher starting amount of CSF.



**Preserved information.** The Volcano plot comparing neat CSF and two different ENRICH-iST preparations (20  $\mu$ L and 200  $\mu$ L) shows the completely depleted and are still detectable and quantifiable.

## **KEY TAKEAWAYS**

- ENRICH-iST: A new pipeline for inhiah-throughput depth, plasma profiling proteome now applicable for cerebrospinal fluid.
- Enhanced Analysis of Low-Abundance **Proteins: I**he technology improves the detection low-abundance proteins, of enabling deeper proteome coverage up to 1.8 fold increase and can also deal with low concentrated samples.
- Drive **Discovery:** Biomarker ENRICH-iST is highly effective for discovering protein new biomarkers of cerebrospinal fluid proteomic studies.
- Data processing using **ProteoScape**<sup>™</sup> (timsDIA-NN real-time algorithm) enables timsTOF MS-runs, analysis resulting in high accurate data for further downstream analysis.

This antibody-free workflow is not species-specific. It can be applied to a variety of mammalian species, which makes it particularly suitable for preclinical studies.

## **More Information**



Conflict of Interest Disclosure: K.Limm, Z.Demianova, K.Hartinger, and N.A. Kulak are mployed by PreOmics GmbH

www.preomics.com