PREOMES

Unlocking the proteomic potential of FFPE tissues with BeatBox[®] and iST: A xylene-free, high-throughput workflow for in-depth proteome analysis

Brianne Nunez¹; Katharina Limm²; Silvia Wuertenberger²,Katrin Hartinger², Jonathan R Krieger³, Nils A. Kulak² ¹PreOmics Inc., Islandia, NY, USA; ²PreOmics GmbH, Martinsried, Germany; ³Bruker Canada Ltd, Milton, ON, Canada

SPOTLIGHT

LC-MS-based proteomics workflow from sample homogenization to data analysis

- Streamlined FFPE sample preparation workflow combining BeatBox and iST technology
- High-throughput processing of up to 96 samples/day to clean peptides
- Optimized xylene-free approach without deparaffinization LC-MS for based proteomics
- Equivalent workflow performance for FFPE samples compared to fresh frozen tissue
- ProteoScape[™] for Coupled to Bruker data accelerated and improved processing

MATERIALS & METHODS

Input:

- Formalin-fixed, paraffin-embedded (FFPE) mouse tissue: 10-µm sections prepared with either xylene-based deparaffinized tissue or full curls
- Fresh frozen mouse tissue (cardiac muscle, kidney and liver; 1-2 mg tissue)

Sample Preparation:

- BeatBox homogenization: BeatBox Tissue Kit 96; 10 min, standard power settings for fresh frozen tissue as "control workflow" or high power settings for FFPE samples)
- Boiling: 1h at 95 °C for FFPE; 10 min at 95 °C for fresh frozen tissue
- iST workflow with optimized washing for FFPE full curls
- workflow: Standard workflow for • iST deparaffinized and fresh frozen tissue; with optimized washing for FFPE full curls

LC-MS analysis & Data processing:

- EASY-nLC[™] 1200 TimsTOF HT
- DIA-PASEF acquisition; 30-min gradient
- ProteoScape[™] Bruker (timsDIA-NN, predicted Library from FASTA)



17-20% for precursors and 9-12% for the respective proteins.



Deep proteome coverage. Protein quantification demonstrated a comparable dynamic development or regulation (based on GO-terms) of heart, kidney or liver are range with approx. 4 orders of magnitude and a similar proteomic depth for full curl shown in the heat map, demonstrating that biological information is preserved FFPE tissue and fresh frozen tissue (three valid values in four replicates).



Preserved information. An evaluation of the biological information conserved between fresh-frozen and FFPE tissue was performed by applying statistical analysis and fitering for proteins significantly regulated for each tissue type. This resulted in 1259 proteins in FF and 1140 proteins in FC, showing a differentially expression between tissue types, but similar protein abundance level within a tissue specimen. A selection of proteins described for the between FFPE and fresh-frozen tissue.





- Optimized solution combining the BeatBox and iST technology provides a simple, fast, and robust way to process FFPE tissue for LC-MS based proteomics.
- Preparing FFPE samples with BeatBox-iST technology saves valuable time by reducing the number of processing steps and eliminating the need for deparaffinization with xylene.
- An innovative FFPE workflow that offers similar overall performance and proteomic depth to fresh frozen tissue is perfectly suited large-scale retrospective for studies with FFPE tissue.
- A more detailed analysis revealed biological conserved information on protein level for selected proteins within a tissue type.

This BeatBox-iST workflow simplifies large-scale retrospective proteomic studies by providing a xylene-free, robust and high-throughput solution for FFPE tissue samples.

More Information & Contact



Contact: Brianne.Nunez@preomics.com

Canada Ltd. Milton. ON. Canada

www.preomics.com