


Overview

OBJECTIVES • Proof-of-concept study for the integration of previously developed untargeted metabolomic and lipidomic workflows, including sample preparation, LC-MS analysis, data processing and identification.


EXPERIMENTAL • Liquid-liquid extraction of serum and urine samples with dichloromethane and methanol; the aqueous layer was used for metabolomics and the organic layer, for lipidomics. Both layers were analyzed by LC-MS, followed by data alignment, filtering, identification, normalization and statistics.

RESULTS • The integrated workflow provided suitable results and reduced the sample volume requirement, as well as overall analysis time.


Introduction



MOTIVATION: (1) Small sample volumes require integration of preparation techniques • (2) Integrated sample preparation allows for high throughput • (3) Overlap between detection by different workflows is reduced.



Methods



- Metabolomics¹**
- MPA: 0.1% formic acid in water
 - MPB: 0.1% formic acid in acetonitrile
 - 22 min gradient, 8 min equilibrium, 30° C
- Lipidomics**
- MPA: 10 mM NH₄COOH in 50:40:10 methanol/acetonitrile/water
 - MPB: 10 mM NH₄COOH in 95:5 isopropanol/water
 - 24 min gradient, 10 min equilibrium, 45° C

Experimental design

Dilution test

Final volume fixed at 30 µL, sample volume varied

Evaluated dilution factors: 1:10, 1:5, 1:2.5 and 1:1

Metabolomics

6 extraction replicates: aqueous fraction

Positive and negative ionization

Lipidomics

6 extraction replicates: organic fraction

Positive and negative ionization

Data processing

Alignment

- Bruker MetaboScape 4.0
- Mass re-calibration, peak picking, alignment, deisotoping
- Filtering (>80% of injections per group)

Normalization

- Lipidomics: internal standardization by 14 deuterated lipids (lipid class and retention time match), summed intensity ratio and auto-scaling
- Metabolomics: summed intensity and auto-scaling




Statistics

- MetaboAnalyst 4.0
- Volcano plots, dendrogram, PCA

Identification





- Metabolomics**
- Positive:** m/z error <5 mDa, mSigma <100, MS/MS score >200
 - Bruker HMDB metabolite library (retention time and MS/MS)¹
 - MoNA² and MSDial LipidBlast³ MS/MS libraries
 - Putative:** MyCompoundID (m/z tolerance of 5 mDa)⁴
- Lipidomics**
- Positive:** m/z error <5 mDa, mSigma <100, MS/MS score >100
 - MSDial LipidBlast, MoNA and HMDB MS/MS libraries^{1,2,3}
 - Putative:** LipidMaps⁵ (m/z match; tolerance of 5 mDa)

Extraction of lipid internal standards



Results

Lipidomics



Conclusions

- The integrated workflow for lipidomics and metabolomics was suitable for the analysis of serum and urine samples.
- The metabolomics methodology allowed de identification of X metabolites.
- The combined sample preparation resulted in the detection of more than 17000 features for serum and 16000 for urine from 12.0 µL of sample.
- The lipidomics methodology allowed the identification of 2925 lipids for serum and 1681 for urine.

Acknowledgements

This work was supported by grants from the Natural Sciences and Engineering Research Council of Canada, Canadian Institutes of Health Research (CIHR), Canada Research Chairs (CRC), Genome Canada, Canada Foundation for Innovation and Alberta Innovates. The authors also acknowledge the help of Xiaohang Wang (University of Alberta, Canada).

References

- <https://www.bruker.com/products/mass-spectrometry-and-separations/ms-software/metabolomics-solution.html>
- Mass Bank of North America, <http://mona.fiehnlab.ucdavis.edu/>
- Kind T, Liu KH, Lee DY, deFelicis B, Meissen JK, Fiehn O, Nature Methods 2013, 10, 755-58.
- LipidMaps database, <https://www.lipidmaps.org/>
- MyCompoundID database, <http://www.mycompoundid.org>