

# Royal jelly's phenolic profile via UPLC-VIP-HESI-TIMS-QTOF-MS: A thorough characterization following a multivariate optimization approach

Elena S. Nastou<sup>1</sup>; Dafni V. Preza-Mayo-Katakis<sup>1</sup>; Panagiotis-Loukas P. Gialouris<sup>1</sup>; Carsten Baessmann<sup>2</sup>; Nikolaos S. Thomaidis<sup>1</sup>

<sup>1</sup>National and Kapodistrian University of Athens, Department of Chemistry, Laboratory of Analytical Chemistry, Athens, Greece; <sup>2</sup>Bruker Daltonics GmbH & Co. KG, Bremen, Germany

## Introduction

A growing scientific interest in the analysis of royal jelly has been developed...

... due to its nutritional and financial significance.



Evidence relates its **positive health impact** to the high content of specific bioactive components, especially **phenolic compounds**<sup>[1]</sup>.

### Incorporation of VIP-HESI & TIMS in LC-HRMS workflows:

Due to the active exhaust

- Higher sensitivity & robustness

**VIP-HESI**

Due to a moving vacuum-insulated probe

- Reduced thermal degradation of sensitive molecules [2]

- Reduced source contamination & memory effects in source

- Addition of a new dimension of separation


**4D- METABOLOMICS**

**TIMS**

- Additional identification confidence in wide-scope target screening [3]

## Methods

### Sample preparation




- Instrumentation: UPLC-VIP-HESI-TIMS-QTOF-MS
- m/z range: 50-1100
- 1/k0 range: 0.40-1.40 Vs cm<sup>-2</sup>
- Ionization mode: Negative
- Acquisition modes: bbCID & PASEF




## VIP-HESI Source Optimization

### Step 1: Sprayer's position and Nebulizer gas pressure

- 9 Sprayer positions inside VIP-HESI



Spiked royal jelly sample with 15 phenolic compounds  
Equivalent results for all!



- 5 Nebulizer gas pressure values

1.5, 2.0, 2.5, 3.0, **3.5 bars**

↑ The pressure of the nebulizer gas ↑ the intensities of all the analytes.

● **Near** positions give **higher peak areas** than respective **Far** positions.

● **Up** positions give **higher peak areas** than respective **Down** positions.

### Step 2: Multivariate optimization of five VIP-HESI parameters

A: Capillary voltage\*

B: Dry gas flow


C: Dry gas temperature

D: Probe gas temperature\*

E: Probe gas flow\*


\*Statistically Significant Parameters

Optimized Parameters:  
Capillary voltage: 2500V  
Dry gas flow: 10 L min<sup>-1</sup>  
Dry gas temperature: 230 °C  
Probe gas flow: 5 L min<sup>-1</sup>  
Probe gas temperature: 450 °C



## References


- [1] Giampieri, F.; Quiles, J. L.; Cianciosi, D.; Forbes-Hernández, T. Y.; Orantes-Bermejo, F. J.; Alvarez-Suarez, J. M.; Battino, M. Bee Products: An Emblematic Example of Underutilized Sources of Bioactive Compounds. *J. Agric. Food Chem.* **2022**, *70* (23), 6833–6848.
- [2] Kourtchev, I.; Szeto, P.; O'Connor, I.; Popoola, O. A. M.; Maenhaut, W.; Wenger, J.; Kalberer, M. Comparison of Heated Electrospray Ionization and Nano-electrospray Ionization Sources Coupled to Ultra-High-Resolution Mass Spectrometry for Analysis of Highly Complex Atmospheric Aerosol Samples. *Anal. Chem.* **2020**, *92* (12), 8396–8403.
- [3] Ridgeway, M. E.; Lubeck, M.; Jordens, J.; Mann, M.; Park, M. A. Trapped Ion Mobility Spectrometry: A Short Review. *International Journal of Mass Spectrometry* **2018**, *425*, 22–35.



HELLENIC REPUBLIC

National and Kapodistrian  
University of Athens

EST. 1837



## Method Validation

- ✓ Linearity
- ✓ Sensitivity
- ✓ Matrix Effect
- ✓ Accuracy
- ✓ Precision


Fit for purpose  
for 55 phenolic  
compounds

	R <sup>2</sup>	LOD (µg g <sup>-1</sup> )	LOQ (µg g <sup>-1</sup> )	ME%	Recovery% (0.25 µg g <sup>-1</sup> )	Recovery% (1.0 µg g <sup>-1</sup> )	Recovery% (5.0 µg g <sup>-1</sup> )
Minimum value	0.98	0.00051	0.0016	-121	80	86	89
Maximum value	0.9992	0.10	0.31	88	101	109	106

	RSD% Repeatability (0.25 µg g <sup>-1</sup> )	RSD% Repeatability (1.0 µg g <sup>-1</sup> )	RSD% Repeatability (5.0 µg g <sup>-1</sup> )	RSD% Intermediate Precision (0.25 µg g <sup>-1</sup> )	RSD% Intermediate Precision (1.0 µg g <sup>-1</sup> )	RSD% Intermediate Precision (5.0 µg g <sup>-1</sup> )
Minimum value	1.0	0.64	1.1	1.3	0.87	1.2
Maximum value	5.9	6.7	5.7	10	12	19

## Target screening



A total of **35 compounds** were identified and quantified in royal jelly samples.

Analyte	C MIN (µg g <sup>-1</sup> )	C MAX (µg g <sup>-1</sup> )	MEAN (µg g <sup>-1</sup> )	MEDIAN (µg g <sup>-1</sup> )	No of samples detected
Apigenin	0.10	1.2	0.34	0.26	22/22
Chrysin	0.11	5.7	1.1	0.62	22/22
Ferulic acid	<LOD	8.6	1.3	0.52	15/22
Galangin	0.21	2.3	0.66	0.61	22/22
Luteolin	0.44	4.7	1.3	1.6	22/22
Naringenin	<LOD	0.79	0.15	0.91	21/22
Pinocembrin	0.15	7.4	1.2	0.38	22/22
Quercetin	0.67	10	2.8	0.68	22/22
Sakuranetin	0.18	14	2.0	0.50	22/22


Mass accuracy ( $\Delta m/z < 2 mDa$ )  
Retention time ( $\Delta RT \pm 0.2 min$ )  
Sigma values (mSigma < 50)  
Qualifier ions  
CCS value ( $\Delta CCS < 2\%$ )

- 18/35 analytes such as Apigenin, Chrysin, Galangin were present in all samples.
- Quercetin had the highest mean concentration. (2.8 µg g<sup>-1</sup>)

## Conclusions

- Seven different parameters were studied and optimized for the first time to facilitate the determination of various phenolic compounds in bee product matrices.
- The optimum parameters of VIP-HESI enhanced the sensitivity of the source, while the addition of TIMS increased the identification confidence.
- Target and suspect screening resulted in the **finger-printing of phenolic compounds of Greek royal jelly**.

## Suspect screening



Database of 130 phenolic compounds

Identification of 85 phenolic compounds



## Future perspectives

This methodology can be extended to assess factors essential to guarantee royal jelly's quality, establish authenticity, and detect fraud cases in samples safe guarding its health promoting properties and allowing beekeepers to expand their businesses obtaining highly competitive prices for their natural product.



Scan me!