


A Novel Approach in Addressing the Challenges of Monitoring Multi-classes of POPs in a Single Run by GC-Ion Mobility-HRMS



Arnd Ingendoh, Carsten Baessmann, Javier Lopez*, Miguel Angel Perez*, Hugo Mueller#, Gauthier Eppe#

Bruker Daltonik GmbH, Bremen, Germany;

*Bruker Espanola, Madrid, Spain;

Mass Spectrometry Laboratory, MoSys Research Unit, Chemistry Department, University of Liège, Belgium

Introduction

Polychlorinated dioxins and furans are bio-accumulative molecules formed during combustion and industrial manufacturing processes. Their analysis is complex due to low regulatory exposure limits and difficult sample matrices. They are persistent organic pollutants (POP), widely found in environmental samples. Severe consequences, even at low exposure concentration, include cancer, reproduction and growth issues, immune system diseases and endocrine effects. Dioxins and furans are mainly quantified in natural matrices by high-resolution sector field MS.

Objectives

1 Demonstrate the capabilities of the GC-APCI-timsTOF for dioxin determination matching the performance criteria of the worldwide standard methods (EPA 1613B, EN 1948) for environmental and food samples.


2 Proof GC-APCI-timsTOF technology as a superior alternative to traditional magnetic sector MS with the benefits of higher analytical flexibility for other POPs and the separation of isomers/isobars for ultra-trace analysis in complex matrices.

3 Develop a broad scope method for POPs screening and quantitation for future regulation compliance.

Analytical performance criteria: EPA 1613B



MS Resolution	>10,000 for HRMS systems
m/z ions	Two exact mass ions to be monitored for each congener
Identification	By retention times and ion ratios. Check RRT native/ ¹³ C
Quantitation	By isotope dilution technique
Calibration	Calibration curve with 5 levels: CS1-CS5
Linearity	Response factor (RF) RSD < 35% coefficient of variation
Ion ratios QC	Tolerance: +/- 15%
Precision	Verification with CS3 level: recoveries RSD
Sensitivity	Minimum level: first calibration point CS1, S/N>10
Carry-over	Analyze a blank after a sample or after CS3




Total of 419 possible congeners. 17 congeners are classified as toxic (due to 2,3,7,8-chlorine substitution).

Dioxins and Furans	Quantitation ion	Confirmation ion	Internal standards	PCBs	Quantitation ion	Confirmation ion	Internal standards
2,3,7,8-TCDF	305.8982	303.9011	[13C12] 2,3,7,8-TCDF	PCB 77	291.8189	289.9218	[13C12] PCB 77
2,3,7,8-TCDD	321.8931	319.8660	[13C12] 2,3,7,8-TCDD	PCB 81	291.8189	289.9218	[13C12] PCB 81
1,2,3,7,8-PeCDF	339.8592	341.8563	[13C12] 1,2,3,7,8-PeCDF	PCB 105	325.8799	327.8771	[13C12] PCB 105
2,3,4,7,8-PeCDF	339.8592	341.8563	[13C12] 2,3,4,7,8-PeCDF	PCB 114	325.8799	327.8771	[13C12] PCB 114
1,2,3,7,8-PeCDD	355.8541	357.8513	[13C12] 1,2,3,7,8-PeCDD	PCB 118	325.8799	327.8771	[13C12] PCB 118
1,2,3,4,7,8-HxCDF	373.8202	375.8173	[13C12] 1,2,3,4,7,8-HxCDF	PCB 123	325.8799	327.8771	[13C12] PCB 123
1,2,3,6,7,8-HxCDF	373.8202	375.8173	[13C12] 1,2,3,6,7,8-HxCDF	PCB 126	325.8799	327.8771	[13C12] PCB 126
1,2,3,7,8,9-HxCDF	373.8202	375.8173	[13C12] 1,2,3,7,8,9-HxCDF	PCB 156	359.8410	361.8381	[13C12] PCB 156
2,3,4,6,7,8-HxCDF	373.8202	375.8173	[13C12] 2,3,4,6,7,8-HxCDF	PCB 157	359.8410	361.8381	[13C12] PCB 157
1,2,3,4,7,8-HxCDD	389.8152	391.8123	[13C12] 1,2,3,4,7,8-HxCDD	PCB 167	359.8410	361.8381	[13C12] PCB 167
1,2,3,6,7,8-HxCDD	389.8152	391.8123	[13C12] 1,2,3,6,7,8-HxCDD	PCB 169	359.8410	361.8381	[13C12] PCB 169
1,2,3,7,8,9-HxCDD	389.8152	391.8123	[13C12] 1,2,3,7,8,9-HxCDD	PCB 189	393.8020	395.7991	[13C12] PCB 189
1,2,3,4,6,7,8-HpCDF	407.8712	409.7784	[13C12] 1,2,3,4,6,7,8-HpCDF				
1,2,3,4,7,8,9-HpCDF	407.8712	409.7784	[13C12] 1,2,3,4,7,8,9-HpCDF				
1,2,3,4,6,7,8-HpCDD	423.7762	425.7733	[13C12] 1,2,3,4,6,7,8-HpCDD				
OCDF	443.7394	441.7423	[13C12] OCDF				
OCDD	459.7343	457.7372	[13C12] OCDD				

List of analyzed compounds: dioxins, furans and dioxin-like PCBs.


Adding the 4th analytical dimension: CCS






Methods

1 μ L sample separated by GC (35 min run, Restek 60 mx0.25 μ m). GC-APCI source coupled to an ion mobility QTOF (timsTOF Pro 2, Bruker) enabling fast and sensitive quantitative analysis of the different classes of dioxins, furans and PCBs in a single GC/HRMS run. Criteria for validation and quantification of compounds: mass accuracy, retention time, isotope pattern matching, MS/MS qualifiers and collision cross sections (CCS) from ion mobility filtering. Kendrick mass defect plots filtered specific compounds containing Cl or Br from the complete GC/MS chromatogram. TargetScreener 4D software (Bruker) used as a post-processing tool for the automated quantification of the POP.


Results



Dioxins and dioxin-like PCB: linearity and RSD.



Analysis of dioxins in palm fatty acid



Conclusions

- A new era adding the 4th dimension for Dioxins and POPs analysis in complex matrices
- Reduction of potential isobaric interferences in complex matrices
- Separation of isomers/isobaric compounds
- More precise quantitation
- Outstanding sensitivity: Low femtograms level in real samples
- Analyze not only Dioxins but a broad scope of POPs in complex matrices
- Prepared for future regulation compliance

timsTOF Pro 2