


Characteristics of MALDI-imaging on a new dual ion source QTOF with TIMS separation



HUPO 2019



Arne Fuetterer¹, Juergen Suetering¹, Janina Oetjen¹, Niels Goedecke¹, Stephanie Kaspar-Schoenefeld¹, Sven Meyer¹, Torsten Mueller², Christian Marsching³, Carsten Hopf³, Jeroen Krijgsveld², Jan H. Kobarg⁴, Nikolas Kessler¹, Wiebke Timm¹, Aiko Barsch¹, Heiko Neuweber¹, Dennis Trede⁴, Shannon Cornett⁵, Alice Ly¹, Lucy Woods¹, Oliver Raether¹, Jens Höhdorf¹

¹ Bruker Daltonik GmbH, Bremen, Germany; ² Deutsches Krebsforschungszentrum (DKFZ), Heidelberg, Germany; ³ Center for Mass Spectrometry and Optical Spectroscopy (CeMOS), Mannheim Technical University, Mannheim, Germany; ⁴ SCiLS Lab, Bremen, Germany; ⁵ Bruker Daltonics Inc., Billerica, MA, USA


Introduction

MALDI Imaging has emerged as a technique with a broad range of applications. We present here the timsTOF fleX, a system consisting of a timsTOF Pro QTOF mounted with a fully integrated high throughput, high spatial resolution MALDI source and stage, enabling measurement of analytes at high mass resolving power with both high mass accuracy and high lateral resolution. The instrument has both an ESI and MALDI source (figure 1), making it ideal for SpatialOMx studies. Additionally, we present a novel software workflow for the identification of signals using data from the high-speed, high spatial resolution timsTOF fleX instrument (figure 2).

Methods



Results



LC-PASEF






Figure 4: SpatialOMx workflow for mouse brain lipids. Mouse brains were dissected, with sections from one hemisphere prepared for MALDI Imaging; lipids were extracted from the other hemisphere (25µg/µl) for LC-MS/MS analysis. Top left panel shows spatial segmentation from SCiLS Lab 2020 of a mouse brain section coated in ZSA matrix [3], and measured in negative mode at 20µm spatial resolution. Key brain regions of interest were exported: cortex (red), corpus callosum (purple), cerebellar white matter and hindbrain (teal) and grey matter (orange). Lower left panel shows 5 replicate LC-PASEF measurements in negative mode of brain lipid extracts from the timsTOF fleX from a 10 µl injection [4]. LC-MS/MS measurement of total brain extract contained 200 compounds. The LC-MS/MS data and MALDI Imaging ROIs were combined in MetaboScape 5.0. Top right panel shows 20 out of 49 annotated features from the combined datasets. Lower right panel shows an example of an annotated lipid in SCiLS Lab 2020.

Conclusions

- The timsTOF fleX delivers uncompromised shotgun proteomics performance making use of PASEF technology while providing MALDI imaging capabilities with 20 µm pitch spatial resolution and 15 pixel/sec speed.
- Addition of the MALDI source to the timsTOF Pro does not compromise proteomics performance
- Lipid annotation automated with MetaboScape 5.0 featuring T-ReX² feature extraction technology and AQ scoring. Visualization of annotated signals with SCiLS Lab completing the intuitive workflow.
- This workflow enables SpatialOMx by providing annotations to the morphological topography

Acknowledgements

We would like to thank Gaby Traue and Eric Mahu from Bruker Daltonik GmbH, and Alexander Geisel, CeMOS, Mannheim, Germany for excellent technical support.

References

- [1] Erich K et al. (2018). Molecular & Cellular Proteomics. doi: mcp.RA118.000980;
- [2] F. Meier et al. MCProt. December 1, 2018, 17 (12) 2534-2545
- [3] Fülöp A et al., Anal. Chem. 2013 Oct 1;85(19):9156-63
- [4] ???

timsTOF fleX