



## ● GC-APCI-QTOF MS for Detection of PCBs and PCDDs: Classic innovation to meet today's trace detection and target quantitation requirements

The reliable detection and quantitation of low-to-medium polarity persistent organic pollutants in food or environmental samples as required by governmental regulations can be readily made via GC-APCI-QTOF MS. High sensitivity is supported by the soft ionization characteristics of the APCI source, preserving molecular ions. Using QTOF MS for the detection and subsequent fragmentation of these ions provides high resolution, high value data for broad screening workflows. The analytical performance of this workflow is demonstrated to meet specific US EPA criteria both in standard solutions and in a complex fish tissue matrix.

### Introduction

The development of analytical instrumentation to detect and confidently identify compounds of interest has progressively

evolved along with the wide array of technological advances over the last decades. Certain "classic" techniques have more recently been recognized as high value elements to modern

analytical workflows. Gas chromatography (GC) has been utilized as a means of sample separation since the late 1950s. Its applicability in analytical schemes is well established,

*Keywords:*  
APCI, GC, QTOF MS, dioxins, persistent organic pollutants (POPs), pesticide detection, food and feed safety, environmental screening




Figure 1: Principle of APCI. Ionization is made via charge transfer when using a "dry" source (top) and via protonation when using a "modified" source (bottom). APCI using a dry source is more suitable for non-polar compounds, while the modified source is more suitable for polar compounds.




Figure 2: GC-APCI-QTOF schematics

benefiting from versatility of column chemistries and precision control of separatory temperature gradient programming, to support efficient and reproducible chromatographic resolution. Following compound separation, the identification and/or quantitation of target compounds typically relies on mass spectrometry (MS), along with some means of effluent ionization, traditionally under vacuum, to permit direct connection between the two instruments.

In recent years, the utility of GC coupled with atmospheric pressure chemical ionization (APCI) has been reintroduced. GC-APCI was developed in the mid-1970s by Dr. Evan Horning and his coworkers at The Institute for Lipid Research of the Baylor College

of Medicine [1,2]. Via a transfer line, GC effluent is introduced into the APCI source chamber, where ionization is initiated (most typically) by a corona discharge. The APCI source chamber (Figure 1) may be used "dry," simply with  $N_2$  carrier gas, forming highly reactive  $N_2^{+\bullet}$  or  $N_4^{+\bullet}$  species which predominantly result in the generation of  $M^{+\bullet}$  molecular ions. The use of a dopant enriched nitrogen favors the formation of protonated molecular ions  $[M+H]^+$ . The dry approach is more suitable to non-polar target compounds, while the doped source can provide improved sensitivity for more polar compounds. Negative ions may also be generated in the APCI source by electron transfer or ion clustering processes.

In order to fully capitalize on the advantages of this separatory and ionization approach, a high resolution, high mass accuracy MS platform is required to anchor the analytical system (Figure 2). Quadrupole time-of-flight (QTOF) mass spectrometry meets this demand, while also enabling the detection of an unlimited number of compounds with high sensitivity across a broad dynamic range. Whether for targeted or discovery workflows, an extremely rich data set can be collected in minutes using alternating full scan and MS/MS acquisition modes.

In many GC-MS analyses, the more commonly used electron (impact) ionization (EI) results in significant target fragmentation within the ion source. Although extensive fragmentation can be valuable in the creation of pattern-based spectral libraries, lower resolution MS systems or co-eluting compounds may restrict confident identification. In contrast, APCI is a lower energy, "soft" ionization technique that maintains the parent molecular ion(s) and limits target fragmentation. This provides several key analytical advantages, including the high sensitivity necessary to detect lower abundant components (as mass accurate  $M^{+\bullet}$  or  $[M+H]^+$  species) within complex samples, and relatively simple fragmentation patterns (Figure 3). APCI-generated ions may then be subjected to controlled fragmentation within the QTOF MS to confirm target identity or characterize novel compounds. APCI is amenable to compounds of low to moderate polarity, up to ~1 kDa, including a much lower polarity range than that of traditional EI (Figure 4). These features support its use for the detection of broad classes of drugs, pesticides, and other environmental pollutants in toxicology, environmental, and food safety screening workflows.




Figure 3: EI-MS and GC-APCI-QTOF MS spectra of endrin ( $C_{12}H_8Cl_6O$ ), a banned herbicide. Using EI (top), the molecule is highly fragmented, such that no molecular ion is detected, and nominal masses are reported. Using GC-APCI-QTOF (bottom), the  $M^+$  ion is prominent, the MS/MS fragmentation pattern (inset, collected simultaneously) is simplified, and all peaks are detected with high mass resolution.

Rule name	Type	Value
1 <input checked="" type="checkbox"/> carryOverLimit	Implicit	300.00
2 <input type="checkbox"/> detectionLimit	Implicit	1.0
3 <input checked="" type="checkbox"/> estimatedLowValue	Implicit	
4 <input checked="" type="checkbox"/> falseNegative	Implicit	
5 <input checked="" type="checkbox"/> lowerLimitOfQuantitation	Implicit	0.50
6 <input checked="" type="checkbox"/> maximumDeviation	Implicit	20.00
7 <input checked="" type="checkbox"/> maximumDeviationQC	Implicit	20.00
8 <input checked="" type="checkbox"/> minimumCorrelation	Implicit	0.99
9 <input checked="" type="checkbox"/> quantityAboveMaxQuantityOfCalibration	Implicit	
10 <input checked="" type="checkbox"/> quantityBelowMinQuantityOfCalibration	Implicit	
11 <input checked="" type="checkbox"/> quantityExceedsReportingLimit	Implicit	1.00
12 <input checked="" type="checkbox"/> responseRange	Implicit	10.00
13 <input checked="" type="checkbox"/> retentionTimeDrift	Implicit	10.00
14 <input checked="" type="checkbox"/> signalBlank	Implicit	0.50
15 <input checked="" type="checkbox"/> standardDeviationResponseFactor	Implicit	35.00
16 <input checked="" type="checkbox"/> upperLimitOfQuantitation	Implicit	1,000.00

Figure 5: TASQ software facilitates the definition of target evaluation criteria according to any regulatory directive requirements. Evaluation criteria for US EPA 1613B are shown.

This note describes the use of GC-APCI-QTOF MS for the determination of dioxin-like polychlorinated biphenyls (PCBs), polychlorinated dibenzodioxins (PCDDs), and polychlorinated dibenzofurans (PCDFs) according to analytical performance criteria of method EPA 1613B [3] in solution and in prepared fish tissue. These dioxin/dioxin-like compounds are considered Persistent Organic Pollutants (POPs), with detrimental and long-lasting effects on human and animal health [4]. Many are highly toxic and are subject to

numerous governmental regulations across the globe. Unambiguous detection, method sensitivity, and quantitation capabilities are required for accurate monitoring, and these demands may be met via GC-APCI-QTOF mass spectrometry.

## Materials and Methods

Targeted congeners of dioxin-like PCBs and PCDDs/PCDFs are shown in Table 1, along with their relative toxicities. Native and  $^{13}C$  labelled targets in n-nonane (Wellington




Figure 4: Ionization ability of different ion sources, indicating typical compound amenability and molecular weight ranges.

Laboratories) were analyzed by GC-APCI-QTOF MS at concentrations ranging from 0.5 - 2000 ng/mL. Fish tissue samples were spiked with a subset of the targeted compounds and prepared according to US EPA guidelines using Pressurized Fluids Extraction (PFE). Instrument conditions are detailed in Table 2.

Data screening and quantitation was made within TASQ (Target Analysis for Screening and Quantitation) software (Bruker Daltonics). All analytical performance criteria described in EPA 1613B were included within a custom method (Figure 5). Automatic mass filtering based on the expected molecular formulae was used ( $\pm 3$  mDa). Within the TASQ software, MRSQ scores were automatically determined for each target compound, including values for mass accuracy, retention time, isotopic pattern fidelity, and the detection of (additional) diagnostic ions (Figure 6).

## Results

### Calibration Characteristics

The calibration curves had excellent linearity, with  $R^2 > 0.998$  for all PCBs and  $R^2 > 0.997$  for all PCDDs/PCDFs. The response factor (RF) RSD was  $<13.7\%$  for the PCBs and  $<13.3\%$  for the PCDDs/PCDFs (Figure 7). The EPA S/N criteria for the lowest calibration point ( $S/N \geq 10$ ) was achieved for all targets (Figure 8). A larger

mass deviation resulted in a lower MRSQ score for one targeted PCDD (1,2,3,6,7,8-HxCDD), although still falling within the wider limits established in the method and within the EPA method requirements.

#### % RSD and Recovery

Precision was calculated by the TASQ software for each isotopically labeled compound within the calibration batch, with each batch including five calibration levels. The RSD was <12.8% for all PCDD/PCDF compounds and <7.3% for all the PCB compounds (Figure 9). Using the median calibration concentration for each target, recoveries ranged from 98-115%. The RRT native/<sup>13</sup>C ratios were from 1.00-1.01, and peak resolution was greater than 18,000 for all targets, easily meeting EPA method requirements (Figure 10).

#### LOD and LOQ

To experimentally calculate the signal-to-noise ratio (S/N), the lowest calibration level for PCDDs/PCDFs and PCBs was diluted 10-fold with n-nonane. As in previous analyses, 1 μL was injected. As required within the referenced EPA criteria, two exact mass ions were detected for each target. As examples, the LOD/LOQ calculation for two target compounds are shown in Figure 11. The LOD for these targets was determined to be 50 fg on-column, with an LOQ of 150 fg on-column.

#### Carry-over evaluation

According to the EPA criteria, the analytical system used must be free from sample carry-over as demonstrated through the injection of an appropriate sample blank following sequential analyses. Immediately following the analysis of the calibration verification QC of PCDDs/PCDFs targets (median concentrations), an n-nonane blank was injected.

Table 1: Studied PCDD/PCDF and dioxin-like PCB congeners and their relative toxicities. Toxic Equivalency Factors (TEF) have been established by the WHO [5] to compare toxicity of dioxins relative to the most toxic TCDD (TEF=1).

Congener	WHO TEF
<b>Dioxins (PCDDs)</b>	
2,3,7,8-TCDD	1
1,2,3,7,8-PeCDD	1
1,2,3,4,7,8-HxCDD	0.1
1,2,3,6,7,8-HxCDD	0.1
1,2,3,7,8,9-HxCDD	0.1
1,2,3,4,6,7,8-HpCDD	0.01
OCDD	0.0003
<b>Furans (PCDFs)</b>	
2,3,4,7,8-PeCDF	0.3
2,3,7,8-TCDF	0.1
1,2,3,4,7,8-HxCDF	0.1
1,2,3,6,7,8-HxCDF	0.1
1,2,3,7,8,9-HxCDF	0.1
2,3,4,6,7,8-HxCDF	0.1
1,2,3,7,8-PeCDF	0.03
1,2,3,4,6,7,8-HpCDF	0.01
1,2,3,4,7,8,9-HpCDF	0.01
OCDF	0.0003

IUPAC N°	Congener	WHO TEF
<b>Dioxin-like PCBs</b>		
126	3,3',4,4',5-PeCB	0.1
169	3,3',4,4',5,5'-HxCB	0.03
81	3,4,4',5-TeCB	0.0003
77	3,3',4,4'-TeCB	0.0001
105	2,3,3',4,4'-PeCB	0.00003
114	2,3,4,4',5-PeCB	0.00003
118	2,3',4,4',5-PeCB	0.00003
123	2',3,4,4',5-PeCB	0.00003
156	2,3,3',4,4',5-HxCB	0.00003
157	2,3,3',4,4',5'-HxCB	0.00003
167	2,3',4,4',5,5'-HxCB	0.00003
189	2,3,3',4,4',5,5'-HpCB	0.00003

Table 2: Instrument conditions

<b>Gas Chromatography</b>	
<b>Instrument</b>	Bruker 436 GC
<b>Carrier gas</b>	Helium, 1.5 mL/min constant flow
<b>Injector</b>	1177 split/spitless, 280°C
<b>Injection volume</b>	1 μL, splitless
<b>Insert</b>	2 mm ID straight liner
<b>Column</b>	BR-Dioxin2, 60 m x 0.25 mm, 0.25 μm film thickness
<b>Column temp. (PCBs)</b>	120°C, hold 0.4 min; 190°C at 50°C/min; 300°C at 2.75°C/min hold 0 min; 320°C at 20°C/min, hold 2.20 min
<b>Column temp. (PCDDs/PCDFs)</b>	120°C, hold 5 min; 250°C at 25°C/min; 300°C at 3°C/min hold 15 min
<b>Mass Spectrometry</b>	
<b>Instrument</b>	Bruker impact II UHR QTOF
<b>Source</b>	GC-APCI
<b>Head temp.</b>	300°C
<b>X-line temp.</b>	280°C
<b>Capillary</b>	4,500 V
<b>Corona</b>	8000 nA
<b>Drying gas</b>	1.5 at 175°C
<b>Nebulization pressure</b>	2.4 bar
<b>MS mode</b>	Full Scan

No peaks were detected, meeting the analytical performance criteria established in the method. False positive identifications are highly improbable due to the method's

automated data processing and multi-faceted detection criteria, and the sensitivity of the system eliminates false negative identifications.




Figure 6: Pictorial, color coded depiction of analysis results within TASQ software. MRSQ scoring depends on key parameters of mass accuracy (M), retention time (R), mSigma\* scores of isotopic fidelity (S), and the presence of qualifier diagnostic ions (Q). This multi-feature scoring eliminates false positive and false negative results.

\*mSigma scores use a scale of 0-1000, where lower values indicate a better fit to theoretical isotopic patterns.

### Fish tissue analysis

Three PCDDs/PCDFs were simultaneously identified and quantitated by GC-APCI-QTOF MS in the fish tissue extract. Using alternate GC column temperature gradient conditions (Table 2), twelve PCBs were likewise simultaneously identified and quantitated. All identification criteria were met according to the referenced EPA guidelines (Figures 12 and 13). The high resolution and mass accuracy of this instrument configuration enables high selectivity in complex sample matrices, such as animal tissue.

## Discussion

The detection and quantitation of the PCDDs/PCDFs and dioxin-like PCBs demonstrated within this study supports the use of GC-APCI-QTOF MS in regulated environmental screening workflows. This technique has also been shown to be effective in the analysis of food packaging contaminants [6,7], other classes of pesticides [8], and drugs of abuse [9]. Although only a small, specific group of toxic pollutants were targeted in this study, other known compounds of interest may be retrospectively

sought within the simultaneously collected MS and MS/MS data pool.

Using Bruker's impact II QTOF MS, there is no practical limit to the number of targets that may be sought and identified. This detection capability, with demonstrated high resolution and high sensitivity, has established QTOF MS as a highly valuable analytical tool for both targeted and untargeted screening approaches. Integrated quality control elements within Bruker's TASQ software avoid false positive and false negatives and demonstrate compliance to




Figure 7: Automated generation of individual PCDD/PCDF calibration curves, demonstrating excellent R<sup>2</sup> values and low RSDs for all compounds



governmental regulatory detection criteria. Although sample separation and introduction have been traditionally made via HPLC and sample ionization via electrospray, well-suited for more polar compounds, the utility of this MS platform to detect and quantify a wider array of compound classes and sample types may be dramatically extended with the use of

a GC-APCI interface. The exchange of ion sources can be made without system venting and with virtually no instrument downtime, supporting rapid turn-around times for result reporting. As shown within this study, target screening requirements in environmental, food safety, occupational health, and toxicology settings can be confidently met via GC-APCI-QTOF MS.

## Further reading:

<https://www.bruker.com/products/mass-spectrometry-and-separations/targetscreener-hr.html>

<https://www.bruker.com/products/mass-spectrometry-and-separations/lc-ms/ion-sources/captivespray-nanobooster/overview.html>

PCDDs/PCDFs									
Calibration Data Point <input type="checkbox"/> Batch Results <input type="checkbox"/> Analysis Results <input checked="" type="checkbox"/>									
MRSQ	Analyte Name	RT Score	RT [min]	m/z Score	$\Delta m/z$ [mDa]	Ions Score	Found.Diag.Ions	mSigma Score	Quantity
1	2,3,4,7,8-PeCDF	++	24.26	++	0.65	++	1	++	1.0 ng/ml
2	2,3,7,8-TCDD	++	20.24	++	0.12	++	1	++	1.1 ng/ml
3	2,3,7,8-TCDF	++	19.76	++	0.52	++	1	++	1.8 ng/ml

PCBs									
Calibration Data Point <input type="checkbox"/> Batch Results <input type="checkbox"/> Analysis Results <input checked="" type="checkbox"/> Chromatograms <input checked="" type="checkbox"/>									
MRSQ	Analyte Name	RT Score	RT [min]	m/z Score	$\Delta m/z$ [mDa]	Ions Score	Found.Diag.Ions	mSigma Score	Quantity
1	PCB 105	++	22.97	++	0.25	++	1	++	1561.1 ng/ml
2	PCB 114	++	22.36	++	0.17	++	1	++	155.5 ng/ml
3	PCB 118	++	21.95	++	0.22	++	1	++	3411.3 ng/ml
4	PCB 123	++	21.73	++	0.24	++	1	++	578.4 ng/ml
5	PCB 126	++	24.23	++	0.06	++	1	++	26.2 ng/ml
6	PCB 156	++	24.52	++	0.20	++	1	++	1527.6 ng/ml
7	PCB 157	++	25.34	++	0.09	++	1	++	447.1 ng/ml
8	PCB 167	++	18.82	++	0.36	++	1	++	10.4 ng/ml
9	PCB 169	++	25.45	++	-0.03	++	1	++	145.4 ng/ml
10	PCB 77	++	21.22	++	0.20	++	1	++	72.7 ng/ml
11	PCB 81	++	20.81	++	0.33	++	1	++	12.8 ng/ml
12	PCB189	++	21.94	++	0.18	++	1	++	7.0 ng/ml

Figure 12: Analysis results for fish tissue extract screening by GC-APCI-QTOF MS. Three PCDDs/PCDFs and twelve PCBs were identified and quantitated within the sample. Color-coded MRSQ scoring indicates all targets are matched within the narrow limits defined in the analysis method.




Figure 13: Example mass filtered ( $\pm 3mDa$ ) EIC chromatogram of targeted analysis of fish tissue extract by GC-APCI-QTOF MS. The dioxin 2,3,7,8-TCDD was identified.

## Conclusion

- GC-APCI-QTOF MS is well suited for the quantitative assessment of targeted PCDDs/PCDFs and dioxin-like PCBs in aqueous and animal tissue samples, exceeding all performance criteria according to EPA-1613B requirements.
- Extensive QC parameters for each target compound, including retention time, mass accuracy, ion ratios, and isotopic pattern fidelity, are defined within the TASQ analysis software, with results easily visualized using color-coded schematics.
- GC-APCI-QTOF MS is a robust and easy-to-use technology for the study of many compound classes, including persistent organic pollutants (POPs). In addition to enabling precise target quantitation, the preservation of molecular ions via APCI together with high resolution QTOF MS analysis supports comprehensive suspect screening and unknown compound characterization in many sample types.



## Learn More



You are looking for further Information?  
Check out the link or scan the QR code for more details.

[www.bruker.com/products/mass-spectrometry-and-separations](http://www.bruker.com/products/mass-spectrometry-and-separations)

### References

- [1] Dzidic I, Carroll DI, Stillwell RN, Horning EC (1976). *Comparison of positive ions formed in nickel-63 and corona discharge ion sources using nitrogen, argon, isobutane, ammonia and nitric oxide as reagents in atmospheric pressure ionization mass spectrometry*. Anal Chem, **48**:1763-1768. DOI: 10.1021/ac50006a035
- [2] Horning EC, Carroll DI, Dzidic I, Haegeler KD, Lin S-N, Oertli CU, Stillwell RN (1977). *Development and Use of Analytical Systems Based on Mass Spectrometry*. Clin Chem, **23**:13-21. DOI: 10.1093/clinchem/23.1.13
- [3] U.S. Environmental Protection Agency Method 1613 Revision B. *Tetra- through Octa-Chlorinated Dioxins and Furans by Isotope Dilution HRGC/HRMS* (October 1994).
- [4] [https://www.who.int/docs/default-source/food-safety/dioxins.pdf?sfvrsn=4bcd5f4d\\_1](https://www.who.int/docs/default-source/food-safety/dioxins.pdf?sfvrsn=4bcd5f4d_1)
- [5] [https://www.who.int/foodsafety/chem/2005\\_WHO\\_TEFs\\_ToxSci\\_2006.pdf](https://www.who.int/foodsafety/chem/2005_WHO_TEFs_ToxSci_2006.pdf)
- [6] Cherta L, Portolés T, Pitarch E, Beltran J, López FJ, Calatayud C, Company B, Hernández F (2015). *Analytical strategy based on the combination of gas chromatography coupled to time-of-flight and hybrid quadrupole time-of-flight mass analyzers for non-target analysis in food packaging*. Food Chem, **188**:301-308. DOI: 10.1016/j.foodchem.2015.04.141.
- [7] Canellasa E, Vera P, Domeño C, Alfaro P, Nerina C (2012). *Atmospheric pressure gas chromatography coupled to quadrupole-time of flight mass spectrometry as a powerful tool for identification of non intentionally added substances in acrylic adhesives used in food packaging materials*. Journal of Chromatography A, **1235**:141-148. DOI: 10.1016/j.chroma.2012.02.039
- [8] Gómez-Pérez ML, Plaza-Bolaños P, Romero-González R, Vidal JLM, Frenich AG (2014). *Evaluation of the Potential of GC-APCI-MS for the Analysis of Pesticide Residues in Fatty Matrices*. J Am Soc Mass Spectrom, **25**: 899-902. DOI: 10.1007/s13361-014-0849-4
- [9] Mesihää S, Ketola RA, Pelander A, Rasanen I, Ojanperä I (2017). *Development of a GC-APCI-QTOFMS Library for New Psychoactive Substances and Comparison to a Commercial ESI Library*. Anal Bioanal Chem, **409**: 2007-2013. DOI: 10.1007/s00216-016-0148-y

For Research Use Only. Not for Use in Clinical Diagnostic Procedures.

● **Bruker Daltonics GmbH & Co. KG**      **Bruker Scientific LLC**

Bremen · Germany  
Phone +49 (0)421-2205-0

Billerica, MA · USA  
Phone +1 (978) 663-3660

[ms.sales.bdal@bruker.com](mailto:ms.sales.bdal@bruker.com) – [www.bruker.com](http://www.bruker.com)