

Detecting PFAS beyond the Current Regulatory Request: a Comprehensive Overview of the Contamination in Water by UPHLC-Ion mobility-HRMS

Igor Nikitin, Arnd Ingendoh, Carsten Baessmann, Andrea Kiehne, Eva-Maria Niehaus, David J. Liwara¹, Pim Leonards¹, Jacob de Boer¹, Sicco Brandsma¹
Bruker Daltonik, Bremen, Germany
¹Department of Environment and Health (E&H), Vrije Universiteit Amsterdam, The Netherlands


Introduction

- PFAS: Huge group of > 4,700 registered compounds plus degradation products and precursors
- Lack of reference standards or spectral libraries, plenty of isomers. Targeted analyses of the complete group is impossible.
- Persistent, bio-accumulative and toxic (PBT). Ubiquitous in environment and organisms.

This study:

- Water sampling at 10 sites in the Netherlands
- SPE pre-concentration
- Non-target, unbiased analysis using UHPLC-HRMS combined with trapped ion mobility (TIMS) with "4D" criteria: mass accuracy, isotope pattern fit, MS/MS, CCS
- Processing with TASQ and MetaboScape
- Subsequent identification by (1) library search, (2) screening data against NIST PFAS suspect list (4,700 entries) or (3) denovo unknown identification
- Kendrick mass defect (KMD) filtering for PFAS compounds from background, based on the fluorine content (repeating CF_2 units).

Workflow



Benefits of TIMS




Fig. 1: The additional dimension of trapped ion mobility (TIMS) separates two co-eluting branched isomers, 4-PFOS and 6-PFOS.





Fig. 2: Improvement of S/N values for chromatographic MS peaks by applying TIMS. The largest difference is visible for low concentrations of main interest.

Sampling sites and SPE pre-concentration



Fig. 3: 10 water sampling sites in NL with different infrastructure (industrial, city, recreation).



Sample pre-concentration by SPE. 25 mL surface water effluent + 50 μ L IS. Blanks (50 μ L IS). 4 mL 0.1% NH_4OH in MeOH for conditioning, 4 mL 0.1% NH_4OH in MeOH for elution + reconstitution.

Non-targeted screening for PFAS by library search and suspect lists




Fig. 4: MetaboScape workflow and overview for the PFAS identification in the effluents using the library search and suspect lists.

Name	FORMULA	INCHI
1. 1H-Perfluoro-3,3-bis(trifluoromethyl)hexane	CBHF13	InChI=1S/CBHF13/C1=2(3(6)13,14)15,7(16,17)18(19,10)9(11,12)18(19,20)21(13)H
2. Bis(perfluorooctyl)phosphonic acid	C20HF39O4P	InChI=1S/C20HF39O4P/C21=5(2,7)25,26(29,30)31(33,34)32(35,40)37(45,46)38(49,50)51(53,54)52(55,56)57(60,61)58(62,63)59(64,65)60(66,67)68(72,73,74)29(79,80)28(85,86)31
3. Chloro-2-propenoic acid 3,3,4,4,5,5,6,6-nonafluorohexyl ester	C9HF19FO2	InChI=1S/C9HF19FO2/C1=4(10)19(20)21=2(6)11,12(7)13,14(8)(15,16)17(18)19(1)1-2H2
4. Bicyclo[2.2.1]hept-2-ene, 5,5,6-trifluoro-6-(heptafluoropropoxy)-	C10HF10O	InChI=1S/C10HF10O/C11=6(12)4=1=2(3)4(7)6(13)21=10(19,20)8(14,15)16(17)18(1)1-2,4-5H,3H2
5. 3,3,4-Trifluoro-4-(heptafluoropropoxy)cyclo[4.1.0]2,5*non-7-ene	C12HF10O	InChI=1S/C12HF10O/C12=13(9)12(13)7=5=2(4)7(9)16(15)10(16)17(11)8(19)20(1)1-2,4-7H,2
6. 5-(Nonafluorobutyl)isocytol[2.2.1]hept-2-ene	C14HF17O2	InChI=1S/C14HF17O2/C1=10(2)4=5=11(10,3)8(2)7(6)9(23)15(10)13(17,18)14(19,20)21(16-7)4,5H
7. Heptafluorobutylcamphor	C14HF17O2	InChI=1S/C14HF17O2/C1=10(2)4=5=11(10,3)8(2)7(6)9(23)15(10)13(17,18)14(19,20)21(16-7)4,5H
8. Perfluorobutylsulfonamide	C4HF9NO2S	InChI=1S/C4HF9NO2S/C5=1(6,3)9,10(11)2(7)8(12)13(14)15(16)17(14,15,16)
9. 1-(3,4-Dichlorophenyl)-2,2,3,3,4,4,4-heptafluorobutanamide	C10H4Cl2F7NO	InChI=1S/C10H4Cl2F7NO/C11=5=2=1(4)3(4)5(12)20=7(21)8(13,14)9(15,16)10(17,18)19(1)1-3H,(H,20,21)
10. 1-(3,5-Dichlorophenyl)-2,2,3,3,4,4,4-heptafluorobutanamide	C10H4Cl2F7NO	InChI=1S/C10H4Cl2F7NO/C11=4=5(12)16(2)=10(2)7(21)8(13,14)9(15,16)10(17,18)19(1)1-3H,(H,20,21)
11. 2-Propenoic acid, 2-methyl-, 2-[[[1-(2-ethylheptafluorooctyl)sulfonyl]heptafluorooctyl]oxy]ethyl ester	C28H28F17NO8S	InChI=1S/C28H28F17NO8S/C1=4-4(5)12(13)16(2)=10(2)7(21)8(13,14)9(15,16)10(17,18)19(1)1-3H,(H,20,21)
12. 2-Propenoic acid, 2-methyl-, 2-[[[1-(2-ethylheptafluorooctyl)sulfonyl]heptafluorooctyl]oxy]ethyl ester	C28H28F17NO8S	InChI=1S/C28H28F17NO8S/C1=4-4(5)12(13)16(2)=10(2)7(21)8(13,14)9(15,16)10(17,18)19(1)1-3H,(H,20,21)
13. 2-Propenoic acid, 2-methyl-, 2-[[[1-(2-ethylheptafluorooctyl)sulfonyl]heptafluorooctyl]oxy]ethyl ester	C28H28F17NO8S	InChI=1S/C28H28F17NO8S/C1=4-4(5)12(13)16(2)=10(2)7(21)8(13,14)9(15,16)10(17,18)19(1)1-3H,(H,20,21)
14. 2-Propenoic acid, 2-methyl-, 2-[[[1-(2-ethylheptafluorooctyl)sulfonyl]heptafluorooctyl]oxy]ethyl ester	C28H28F17NO8S	InChI=1S/C28H28F17NO8S/C1=4-4(5)12(13)16(2)=10(2)7(21)8(13,14)9(15,16)10(17,18)19(1)1-3H,(H,20,21)

Fig. 5: Detail of the NIST PFAS suspect list containing just very simple information about the name, elemental composition and InChI structure of 4700 PFAS. <https://comptox.epa.gov/dashboard/>

Filtering by KMD plots




Fig. 6: Plotting the Kendrick Mass Defect (KMD) with CF_2 as the base unit versus m/z . Starting from a total of 15,700 detected features, 1181 potential PFAS have been filtered (92% data reduction).

Identified PFAS by various methods




Fig. 7: Venn diagram as a summary of the identified PFAS by the various identification steps. From 15,700 extracted features in the LC/MS/MS run, the different ID steps lead to:

- ID step 1: 23 PFAS could be annotated by library searches in the Bruker PFAS library (green) and NIST PFAS library (blue) with an overlap of 6 compounds.
- ID step 2: Additional 33 PFAS were identified by the comparison with the NIST PFAS suspect list (purple).
- ID step 3: KMD was used to filter the remaining features systematically for potential PFAS. 2 more PFAS were found via the denovo identification workflow (orange).

In total, 58 PFAS were found and identified by this holistic approach, providing an overview of the total content of PFAS in each sample. All steps represent tentative results, the post-processing is still work-in-progress, and more PFAS are expected to be found after a complete evaluation of the data from all sampling sites.





Fig. 8: Various PFAS species with different intensities at the sampling sites. Dordrecht sticks out for PFOA, e.g., while Tilburg has a higher MeFhSAA concentration. As well, new PFAS precursors have been found which are not screened yet in legal directives.



Criterion	Score
MS/MS	814
$\Delta m/z$	-0.37
ΔCCS	5.8
mSigma	6.5
rt	8.96
CCS	166.3

Fig. 9: Example for the denovo ID step 3. PFAS H2-U-PFOS was identified by the MetaboScape tools SmartFormula, Compound Crawler, in-silico fragmentation (MetFrag) and CCS prediction. This unsaturated PFAS has already been identified before in the environment, but a standard is not commercially available yet. known

Conclusion

- The 4D scoring of exact mass, isotope pattern match, MS/MS data and CCS values provides a reliable basis for confident identification.
- Different workflows including (1) library search, (2) comparison with suspect lists and (3) denovo identification have been successfully applied.
- In total, 58 PFAS were found in the water samples, with distinct differences for the various sampling sites which point to individual environmental circumstances and origins for PFAS.
- This wide-scope screening of real-life samples proved to be a comprehensive approach for a fast and efficient PFAS identification and will highly assist in the understanding of the chemical universe of PFAS in the environment and for protecting the environment, wildlife, and human health.

timsTOF Pro 2